



easy*PGX*[®]

Ready to yo**use**

diatech
pharmacogenetics





READY TO USE

Reagents delivered in 8-well strips preloaded with a complete master mix




EASY TO USE

No need for freezing, thawing or pipetting on ice and the few remaining pipetting steps minimize the risk of errors or contamination



HIGH SENSITIVITY

Limit of detection as low as 0.5%



FLEXIBLE SAMPLE REQUIREMENT

Low DNA or RNA input from a variety of sources, including FFPE and plasma



TURNAROUND TIME

From tissue to result in less than 3 hours with only 10 minutes of hands-on time




FLEXIBLE USE

Possibility to run multiple targets all in one experiment: same thermal profile for the somatic assays and color coded reagents



AUTOMATIC DATA ANALYSIS

Import data into Easy PGX dedicated software and get results



TRANSPORT AND STORAGE AT ROOM TEMPERATURE

Complete master mix in a dry format, stable at room temperature



STANDARDS INCLUDED

Positive and negative controls included for validation of each experimental session



AUTOMATIC EVALUATION

Inhibitors or pipetting errors detected by the internal reaction control



REGULATORY

Kits have been designed, developed and validated in accordance with the Directive 98/79/EC on in vitro diagnostic medical devices



QUALITY ASSURANCE

Manufactured under ISO 13485

Easy_{PGX}[®] ready KRAS
cat. no. RT021 (48 test, CE IVD)

MAIN FEATURES

Detection of the main mutations of exon 2 (codons 12, 13), of exon 3 (codons 59, 61) and of exon 4 (codons 117, 146) of the KRAS gene.

Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready BRAF
cat. no. RT022 (48 test, CE IVD)

MAIN FEATURES

Detection of the main mutations of codon 600 of the BRAF gene.

Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready EGFR
cat. no. RT023 (48 test, CE IVD)

MAIN FEATURES

Detection of the main mutations of exons 18, 19, 20, 21 of the EGFR gene.

Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready NRAS
cat. no. RT024 (48 test, CE IVD)

MAIN FEATURES





Detection of the main mutations of exon 2 (codons 12, 13), of exon 3 (codons 59, 61) and of exon 4 (codons 117, 146) of the NRAS gene.

Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and plasma*.

Easy PGX system: from tissue to result in less than 3 hours

| Workflow | One step tissue lysis | PCR setup | PCR run | Data analysis |
|---------------|---|---|---|---|
| Total time |  60 min Add the Easy _{PGX} Dep solution and extraction reagents to the sample |  < 5 min Add extracted samples to ready-to-use 8-well strip |  90 min Load the 8-well strip onto the thermal-cycler and start the run |  < 1 min Import raw data into the Easy _{PGX} Analysis software |
| Hands-on time | < 5 min | < 5 min | < 1 min | < 1 min |

Easy_{PGX}® ready ALK, ROS1, RET, MET
cat. no. RT025 (48 test, CE IVD)

MAIN FEATURES

Detection of the main chromosomal translocations involving ALK, ROS1, RET and the MET exon 14 skipping.

Each mix allows the co-amplification of one or more fusions plus an endogenous control gene.

STARTING MATERIAL

RNA from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues and cytological samples.

Easy_{PGX}® ready DPYD
cat. no. RT026 (48 test, CE IVD)

MAIN FEATURES

Detection, by allelic discrimination, of the DPYD gene polymorphisms DPYD*2A (IVS14+1G>A, c.1905+1G>A), DPYD*13 (c.1679T>G), DPYD D949V (c.2846A>T) and DPYD IVS10 (c.1129-5923C>G), associated with the toxicity due to the treatment with Fluoropyrimidines.

Each mix allows the co-amplification of the mutant sequence as well as the wild-type sequence.

STARTING MATERIAL

DNA from whole blood.

Easy_{PGX}® ready UGT1A1
cat. no. RT027 (48 test, CE IVD)

MAIN FEATURES

Detection, by allelic discrimination, of the UGT1A1 gene polymorphisms UGT1A1*36 (TA)5, UGT1A1*1 (TA)6, UGT1A1*28 (TA)7 and UGT1A1*37 (TA)8, associated with the toxicity due to the treatment with Irinotecan.

Each mix allows the co-amplification of the target polymorphisms plus an endogenous control gene.

STARTING MATERIAL

DNA from whole blood.

Easy_{PGX}® ready THYROID
cat. no. RT028 (48 test, CE IVD)

MAIN FEATURES

Detection of the main mutations of exon 2 (codons 12,13), of exon 3 (codon 61) of the KRAS, NRAS, HRAS genes and of the codons 600 and 601 of the BRAF gene.

Each mix allows the co-amplification of one or more mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed paraffin-embedded (FFPE) tissues, and cytological samples.

Color coded: multiple targets in one run



Easy_{PGX}[®] ready FL-DNA
cat. no. RT029 (48 test, CE IVD)

MAIN FEATURES

Absolute quantitative detection of Fluorescence
Long DNA from faecal specimens.
Each mix allows the co-amplification of the target
DNA (APC and TP53) plus an exogenous control gene.

STARTING MATERIAL

DNA from faecal specimens.

Easy_{PGX}[®] ready EGFR PLUS
cat. no. RT030 (48 test, CE IVD)

MAIN FEATURES

Detection of T790M and C797S (c.2389 T>A, c.2390 G>C) of the EGFR gene.
Each mix allows the co-amplification of one or more
mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed
paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready IDH 1-2
cat. no. RT031 (48 test, CE IVD)

MAIN FEATURES

Detection of the main mutations of IDH1 gene
(codons 105 and 132) and IDH2 gene (codons 140 and
172).
Each mix allows the co-amplification of one or more
mutated alleles plus an endogenous control gene.

STARTING MATERIAL

DNA from fresh, frozen, formalin-fixed
paraffin-embedded (FFPE) tissues and plasma*.

Easy_{PGX}[®] ready THYROID FUSION
cat. no. RT032 (48 test, CE IVD)


MAIN FEATURES

Detection of the chromosomal translocations
involving RET/PTC1: CCDC6-RET; RET/PTC2:
PRKAR1A-RET; RET/PTC3: NCOA4-RET and
PAX8/PPARG.
Each mix allows the co-amplification of one or more
fusions plus an endogenous control gene.

STARTING MATERIAL

RNA from fresh, frozen, formalin-fixed paraffin-
embedded (FFPE) tissues and cytological samples.

Real time reaction monitoring



* Please note that extraction from plasma is sold separately (cat.n. H8040)

EasyPGX® ready MSI cat. no. RT033 (48 test, CE IVD)

MAIN FEATURES

Detection of 8 mononucleotide "quasi - monomorphic" markers: BAT-25, BAT-26, NR-21, NR-22, NR-24, NR-27, CAT-25 and MONO-27 by Real Time PCR and subsequent analysis of the targets based on the denaturation profile. The test allows, accurately and with reduced "hands-on time", to detect the microsatellite instability in tumor samples.

STARTING MATERIAL

DNA from fresh, frozen, formalin fixed paraffin embedded (FFPE) tissue. Comparison with normal tissue or blood is not necessary for the analysis of results.

Helix® circulating Nucleic Acid cat. no. H8040 (50 test, CE IVD)

MAIN FEATURES

The kit allows the manual extraction of circulating free DNA (cfDNA) from plasma. The kit Helix® Circulating Nucleic Acid, in association with the kit EasyPGX® ready EGFR, enables the mutational analysis of EGFR gene in the circulating tumour DNA (liquid biopsy) when the tumour tissue is not evaluable, according to the EMA/129677/2014 recommendations of September 25th 2014.

DNA capture by silica membrane and vacuum-based system. The system to concentrate the final eluate up to 3 times is included in the kit.

STARTING MATERIAL

1-5 ml of fresh or frozen plasma.

TURN AROUND TIME

3 hours.

Automatic data analysis software

File Options Tools

Experiment
Name: Administrator
Experiment description:
Add description here...

Platform and application
Select instrument:
RT800-96_EasyPGX qPCR instrument 96
Select kit:
RT023_EasyPGX ready EGFR

Import data and analysis
Select data:
CARun120170608_RT023_EGFR.xlsx
[Open] [Start] [Save report] [Reset]

Analysis of the reaction controls:

| NL | Name | Dye | G210 | | | T790M | | | S76H | | | rs20136 | | | L858R | | | T861Q | | | rs4546 | | | CTSL | | |
|----|----------|-----|------|-----|--------|-------|-----|--------|------|------|--------|---------|------|--------|-------|------|--------|-------|-----|--------|--------|------|--------|------|------|--------|
| | | | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result |
| 1 | Positive | FAM | 24.4 | 673 | OK | 19.4 | 934 | OK | 20.2 | 2236 | OK | 21.8 | 1043 | OK | 20.6 | 2302 | OK | 21.0 | 830 | OK | 21.0 | 1520 | OK | 24.1 | 1142 | OK |
| 2 | Negative | FAM | 6 | OK | | 2 | OK | | 6 | OK | | 46 | OK | | 5 | OK | | -1 | OK | | 34.2 | 305 | OK | 24.7 | 527 | OK |
| | | HEX | 24.8 | 495 | OK | 24.2 | 698 | OK | 23.4 | 1184 | OK | 23.9 | 584 | OK | 23.3 | 1542 | OK | 24.1 | 749 | OK | 23.4 | 1072 | OK | 24.7 | 527 | OK |
| | | HEX | 7 | OK | | 4 | OK | | 13 | OK | | 9 | OK | | 19 | OK | | 11 | OK | | 10 | OK | | 3 | OK | |
















Analysis of the sample control mix and mutation assay:

| NL | Name | Dye | G210 | | | T790M | | | S76H | | | rs20136 | | | L858R | | | T861Q | | | rs4546 | | | CTSL | | |
|----|---------------------|-----|------|-----|--------|-------|-----|--------|------|------|--------|---------|-----|--------|-------|------|--------|-------|------|--------|--------|------|--------|------|------|--------|
| | | | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result | Cq | AR | Result |
| 1 | EGFR_FFPE_sample 1 | FAM | 4 | WT | | 3 | WT | | 10 | WT | | 16 | WT | | 3 | WT | | 5 | WT | | 3.1 | 314 | WT | 30.7 | 528 | OK |
| | | HEX | 2.3 | 245 | OK | 1.7 | 393 | OK | 1.5 | 519 | OK | 1.2 | 308 | OK | 1.3 | 684 | OK | 1.5 | 551 | OK | 0.9 | 584 | OK | 30.4 | 288 | OK |
| 2 | EGFR_FFPE_sample 2 | FAM | 4 | WT | | 12 | WT | | 2 | WT | | 16 | WT | | 9 | WT | | 0.4 | 554 | MUT | 3.9 | 250 | WT | 30.4 | 882 | OK |
| | | HEX | -0.1 | 387 | OK | -0.4 | 508 | OK | -0.9 | 808 | OK | -1.0 | 301 | OK | -0.8 | 922 | OK | -0.3 | 574 | OK | -0.2 | 641 | OK | 31.3 | 379 | OK |
| 3 | EGFR_FFPE_sample 3 | FAM | 5 | WT | | 1.5 | 695 | MUT | -3 | WT | | 12 | WT | | 7 | WT | | -0.6 | 1528 | MUT | 27.9 | 1113 | OK | | | |
| | | HEX | 0.6 | 412 | OK | -0.1 | 645 | OK | -0.3 | 1060 | OK | -0.6 | 453 | OK | -1.0 | 1412 | OK | -0.5 | 708 | OK | -0.6 | 976 | OK | 28.0 | 545 | OK |
| 4 | EGFR_FFPE_sample 4 | FAM | 3.5 | 525 | MUT | -7 | WT | | -0.2 | 2272 | MUT | 18 | WT | | 9 | WT | | 6.3 | 225 | WT | 27.1 | 1220 | OK | | | |
| | | HEX | 0.6 | 467 | OK | 0.1 | 729 | OK | -0.3 | 1223 | OK | -0.3 | 586 | OK | -0.6 | 1515 | OK | -0.1 | 756 | OK | -0.6 | 1104 | OK | 26.2 | 652 | OK |
| 5 | EGFR_FFPE_sample 5 | FAM | 19 | WT | | 2 | WT | | -13 | WT | | 16 | WT | | 11 | WT | | 6 | WT | | 3.2 | 294 | WT | 30.3 | 861 | OK |
| | | HEX | 0.1 | 418 | OK | 0.2 | 590 | OK | -0.8 | 1150 | OK | -0.7 | 375 | OK | -1.0 | 1231 | OK | -0.5 | 744 | OK | 0.0 | 809 | OK | 29.9 | 476 | OK |
| 6 | EGFR_FFPE_sample 6 | FAM | 7 | WT | | 2 | WT | | 21 | WT | | 16 | WT | | 11 | WT | | 2 | WT | | 5.1 | 284 | WT | 28.9 | 1041 | OK |
| | | HEX | 0.5 | 422 | OK | -0.1 | 651 | OK | -1.3 | 1205 | OK | -0.9 | 500 | OK | -0.2 | 1300 | OK | -0.5 | 752 | OK | -0.7 | 895 | OK | 28.8 | 548 | OK |
| 7 | EGFR_FFPE_sample 7 | FAM | 9 | WT | | 9 | WT | | 13 | WT | | 16 | WT | | 2.3 | 970 | MUT | 8 | WT | | 4.9 | 255 | WT | 28.9 | 886 | OK |
| | | HEX | -0.3 | 552 | OK | -0.5 | 639 | OK | -1.2 | 1094 | OK | -0.9 | 484 | OK | -1.3 | 1372 | OK | -1.3 | 772 | OK | -1.1 | 993 | OK | 28.9 | 445 | OK |
| 8 | EGFR_FFPE_sample 8 | FAM | 4 | WT | | 4 | WT | | 13 | WT | | 14 | WT | | 6 | WT | | -1.4 | 1375 | MUT | 30.1 | 927 | OK | | | |
| | | HEX | 0.4 | 434 | OK | 0.2 | 640 | OK | -0.7 | 1104 | OK | -0.2 | 373 | OK | -0.8 | 1295 | OK | -0.2 | 726 | OK | -0.2 | 903 | OK | 29.3 | 540 | OK |
| 9 | EGFR_FFPE_sample 9 | FAM | 13 | WT | | 6 | WT | | 15 | WT | | 21 | WT | | 29 | WT | | 6 | WT | | 0.0 | 1304 | MUT | 29.6 | 993 | OK |
| | | HEX | 0.6 | 453 | OK | 0.3 | 620 | OK | -0.9 | 586 | OK | -0.1 | 394 | OK | -0.9 | 1329 | OK | 0.0 | 735 | OK | -0.4 | 936 | OK | 28.9 | 585 | OK |
| 10 | EGFR_FFPE_sample 10 | FAM | 4 | WT | | 24 | WT | | 17 | WT | | 14 | WT | | 11 | WT | | 5 | WT | | 4.2 | 290 | WT | 29.3 | 936 | OK |
| | | HEX | 0.3 | 374 | OK | 0.1 | 573 | OK | -1.2 | 1083 | OK | -0.9 | 472 | OK | -1.2 | 1363 | OK | -0.7 | 677 | OK | -1.0 | 987 | OK | 28.4 | 498 | OK |

Warnings

| Code | Description |
|------|---|
| e03 | Possible error in the set up of the reaction / run. It is not possible to analyze the samples (see Troubleshooting) |
| e04 | Possible contamination: it is not possible to analyze the samples (see Troubleshooting) |
| e05 | Suboptimal amount of starting DNA or PCR inhibition (see Troubleshooting) |
| e06 | Excess of DNA. Sample must be diluted with water so that Cq falls in the range indicated (see Troubleshooting) |
| e07 | Not sufficient template / PCR inhibition / mistake during samples dispensation (see Troubleshooting) |
| e08 | Probable excess of DNA. Proceed with the analysis of the sample. |

System ordering information

| Catalog number | Product description | Picture |
|--|---|--|
| RT800-96   | Easy _{PGX} qPCR instrument 96 |  |
| RT800-SW   | Easy _{PGX} analysis software |  |
| RT801   | Easy _{PGX} dry block |  |
| RT802   | Easy _{PGX} centrifuge/vortex 8-well strips |  |
| RT803   | Easy _{PGX} centrifuge/vortex 1.5 ml |  |

Kit ordering information

| Catalog number | Product description | Color coded |
|---|---|---|
| RT021   | Easy _{PGX} ready KRAS (48 test) | Green  |
| RT022   | Easy _{PGX} ready BRAF (48 test) | Purple  |
| RT023   | Easy _{PGX} ready EGFR (48 test) | White  |
| RT024   | Easy _{PGX} ready NRAS (48 test) | Blue  |
| RT025   | Easy _{PGX} ready ALK, ROS1, RET, MET (48 test) | Purple / Blue   |
| RT026   | Easy _{PGX} ready DPYD (48 test) | Purple  |
| RT027   | Easy _{PGX} ready UGT1A1 (48 test) | Clear  |
| RT028   | Easy _{PGX} ready THYROID (48 test) | Black  |
| RT029   | Easy _{PGX} ready FL-DNA (48 test) | White  |
| RT030   | Easy _{PGX} ready EGFR Plus (48 test) | Red  |
| RT031   | Easy _{PGX} ready IDH1-2 (48 test) | Clear  |
| RT032   | Easy _{PGX} ready THYROID Fusion (48 test) | White  |
| RT033   | Easy _{PGX} ready MSI (48 test) | Green / Black   |

For information please contact:

diatech pharmacogenetics

Diatech Pharmacogenetics srl

Via Ignazio Silone 1b - 60035 Jesi (An) Italy

Phone +39 0731 213 243

export@diatechpgx.com

www.diatechpharmacogenetics.com